




TALLER 3



Practical: Sequence your genome, remap and identify variants

Please note that these are example commands and that you might be forced to change them according to your files and paths.

commands are given in “Courier” font

The “\” sign should not be copied, it is here to show that the line continues on the next line...

#login to the main server

```
ssh userxxxx@168.176.61.249
```

#set path for samtools

```
echo export PATH=/data/samtools/bin:$PATH >>.profile
```

#exit session ssh

```
exit
```

#again login to the main server

```
ssh userxxxx@168.176.61.249
```

BWA (mapping)

```
mkdir bwa
```

```
cd bwa
```

#link the reference genome and the reads locally

```
ln -s /data/refgenome/refGenome.fasta refGenome.fasta
```

```
ln -s /data/filtered_reads filtered_reads
```

#index reference for bwa

```
bwa index refGenome.fasta
```

#map reads onto reference



```
zcat filtered_reads/barcode09_filtered.fastq.gz | bwa mem -t 16  
refGenome.fasta - > barcode09.sam
```

BOWTIE (mapping)

```
cd  
mkdir bowtie  
cd bowtie
```

#link the reference genome and the reads locally

```
ln -s /data/refgenome/refGenome.fasta refGenome.fasta  
ln -s /data/filtered_reads filtered_reads
```

build the index of the reference genome

```
bowtie2-build -f refGenome.fasta refgenome
```

map the reads onto the reference save a sam file

```
zcat filtered_reads/barcode09_filtered.fastq.gz | bowtie2 -q -p 16 -  
x refgenome -U - -S barcode09.sam 2> bowtie.log
```

BOTH

#index reference for samtools

```
samtools faidx refGenome.fasta
```

#convert sam to bam

```
samtools view -b -@ 4 -t refGenome.fasta.fai barcode09.sam >  
barcode09_unsorted.bam
```

#sort bam

```
samtools sort -@ 4 -o barcode09_sorted.bam barcode09_unsorted.bam
```

#index bam

```
samtools index barcode09_sorted.bam
```

Extract the consensus

#call variants

```
bcftools mpileup -Ou -f refGenome.fasta barcode09_sorted.bam |  
bcftools call -mv -Oz -o calls.vcf.gz
```

```
bcftools index calls.vcf.gz
```



#normalize indels

```
bcftools norm -f refGenome.fasta calls.vcf.gz -Ob -o calls.norm.bcf
```

#filter adjacent indels within 5bp

```
bcftools filter --IndelGap 5 calls.norm.bcf -Ob -o calls.norm.flt-indels.bcf
```

apply variants to create consensus sequence

```
cat refGenome.fasta | bcftools consensus calls.vcf.gz > consensus.fa
```

Once the mapping is complete, you can obtain some metrics about the alignment as follows:

BWA doesn't provide metrics of the alignment. So, use the next command for counting the number of unmapped reads from the SAM output.

```
awk '{ if ( $3~/\*/ ) {print} }' barcode09.sam | wc -l
```

Bowtie provides some metrics during the alignment process. Look at the entire log file (is very short).

```
more bowtie.log
```

Then, answer the following questions for both single / paired alignment:

- 1. How many reads were mapped onto the reference?***
- 2. How many reads were not mapped onto the reference?***

View with the terminal (? for help, q for quit)

```
samtools tview barcode09_sorted.bam refGenome.fasta
```

Is the reference uniformly covered? Do you observe any region of high coverage? Could you explain why?

Visualize with IGV or Tablet

#download IGV or use the java webstart application
<https://www.broadinstitute.org/software/igv/download>